



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL18 (ATL18)

Product Code	CSB-CF871219DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SZL4
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	KKNGGDADAHDDDGYNLVGVMFGDKEKEEEEICCPICLVEFEAEDAVTHLPRC AHLFHINC IEPWLLRGHLTCPLCRSFVLAPTPTQNVNNAHSSSTLYLSIFFFFCIFLHLLGY L
Source	in vitro E.coli expression system
Target Names	ATL18
Protein Names	Recommended name: RING-H2 finger protein ATL18
Expression Region	30-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.