



Recombinant Escherichia coli O127:H6 Fumarate reductase subunit D (frdD)

Product Code	CSB-CF486511EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7UPX3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Sequence	<pre>MINPNPKRSDEPVFWGLFGAGGMWSAIIAPVMILLVGILLPLGLFPGDALSYER VLAFAQ SFIGRVFLFLMIVLPLWCGLHRMHHAMHDLKIHVPAGKWWFYGLAAILTVVTLIG VVTI</pre>
Source	in vitro E.coli expression system
Target Names	frdD
Protein Names	Recommended name: Fumarate reductase subunit D Alternative name(s): Fumarate reductase 13 kDa hydrophobic protein
Expression Region	1-119
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.