



Recombinant Escherichia coli Fumarate reductase subunit D (frdD)

Product Code	CSB-CF539341ENP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1ITP5
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Sequence	<pre>MINPNPKRSDEPVFWGLFGAGGMWSAIIAPVMILLVGILLPLGLFPGDALSYER VLAFAQ SFIGRVFLFLMIVLPLWCGLHRMHAMHDLKIHVPAGKWWFYGLAAILTVVTLIG VVTI</pre>
Source	in vitro E.coli expression system
Target Names	frdD
Protein Names	Recommended name: Fumarate reductase subunit D Alternative name(s): Fumarate reductase 13 kDa hydrophobic protein
Expression Region	1-119
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.