



Recombinant Inner membrane protein yidG (yidG)

Product Code	CSB-CF364997EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ADL7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MPDSRKARRIADPGLQPERTSLAWFRTMLGYGALMALAIKHNWHQAGMLFWI SIGILAIV ALILWHYTRNRNLMDVTNSDFSQFHVVRDKFLISLAVLSLAILFAVTHIHQLIVFI ERVA
Source	in vitro E.coli expression system
Target Names	yidG
Protein Names	Recommended name: Inner membrane protein yidG
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.