



Recombinant Danio rerio UPF0444 transmembrane protein C12orf23 homolog (si:dkey-227k11.2, zgc:64098)

Product Code	CSB-CF749232DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7T352
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSETNQQRDDLPAFLNDEPPADANKDHPQPQPGMFSRVAGGLFSVTKGAVG ATVGGVAWV GGKSFDLTKTAVTSVASVPAMGVGLVKGGVSAVAGGVAVGSAAVSKVPIGG GKKKDKSD
Source	in vitro E.coli expression system
Target Names	tmem263
Protein Names	Recommended name: UPF0444 transmembrane protein C12orf23 homolog
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.