



Recombinant Human UPF0414 transmembrane protein C20orf30 (C20orf30)

Product Code	CSB-CF853380HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96A57
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMPSTRNLATGIPSSKVKYSRLSSTDDGYIDLQFKKTPPKIPYKAIALATVLFLLIG AFLI IIGSLLLSGYISKGGADRAVPVLIIGILVFLPGFYHLRIAYYASKGYRGYSYDDIPD FDD
Source	in vitro E.coli expression system
Target Names	TMEM230
Protein Names	Recommended name: UPF0414 transmembrane protein C20orf30
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.