



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YGL088W (YGL088W)

Product Code	CSB-CF345262SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53151
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MHEVTRTYFFLFFFLSYKRQINAAFIALDFDPLLFIYFPFLILVLFYNSNANLTAI RNT YSISRLNPSGAFLTHEECGLVLQYIYYWLGLENKFIDLGCNSLSVVCFLADLR VYLRVP G
Source	in vitro E.coli expression system
Target Names	YGL088W
Protein Names	Recommended name: Putative uncharacterized protein YGL088W
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.