



# Recombinant Pig Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial (SDHD)

<b>Product Code</b>	CSB-CF020909PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A5GZW8
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Sequence</b>	DRHTPGWCGVQHIHLSPSHQASSKAASLHWTGERVVSLLLLGLLPAAAYLNPC SAMDYSLA AALTLHGHWGIGQVVTDYVRGDALQKVAKAGLLALSFTFAGLCYFNYHDVGI CKAVAML WKL
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SDHD
<b>Protein Names</b>	Recommended name: Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial Short name= CybS Alternative name(s): CII-4 QPs3 Succinate dehydrogenase complex subunit D Succinate-ubiquinone oxidoreductase cyto
<b>Expression Region</b>	37-159
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.