



Recombinant Pongo pygmaeus NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 11, mitochondrial (NDUFB11)

Product Code	CSB-CF603888EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQJ3
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	ESSFSRTVVAPSAVARKRLPEPTTQWQEDLDPEDENLYEKNPDSHGYDKDPV LDVWNMRL VFFFGVSIILVLGSTFVAYLPDYRMKEWSRREAERLVKYREANGLPIMESNCFD PSKIQL PEDE
Source	in vitro E.coli expression system
Target Names	NDUFB11
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 11, mitochondrial Alternative name(s): Complex I-ESSS Short name= CI-ESSS NADH-ubiquinone oxidoreductase ESSS subunit
Expression Region	30-153
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.