



Recombinant Rickettsia bellii Succinate dehydrogenase hydrophobic membrane anchor subunit (sdhD)

Product Code	CSB-CF631186RAAI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1RHB6
Product Type	Transmembrane Protein
Immunogen Species	Rickettsia bellii (strain RML369-C)
Sequence	MTYDFRAEIVKAKNTGSAKSGSHHWLLQRITAILVLCSLWLLYFTLANKNSDVN IIIWE LKRPNLIPLLIHAVITSLYHAMLGMQVVIEDYISCNKLRNTLIIAVKLFSLTIVAFIVA VFYRG
Source	in vitro E.coli expression system
Target Names	sdhD
Protein Names	Recommended name: Succinate dehydrogenase hydrophobic membrane anchor subunit
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.