



Recombinant Rickettsia bellii Succinate dehydrogenase cytochrome b556 subunit (sdhC)

Product Code	CSB-CF636192RAAI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1RHB5
Product Type	Transmembrane Protein
Immunogen Species	Rickettsia bellii (strain RML369-C)
Sequence	MTKTKQEIYNKRPTSPHLTIYKPKISSTLSILHRMTGVALFFAVSILAWWFILSKF DSNY IKLANCCCIKICLILTSFAWFYHLCNGIRHLFWDIGLGFSSIKAVNLTGWSVVICSV LFT ILLWV
Source	in vitro E.coli expression system
Target Names	sdhC
Protein Names	Recommended name: Succinate dehydrogenase cytochrome b556 subunit Short name= Cytochrome b-556
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.