



Recombinant Drosophila melanogaster UPF0389 protein CG9231 (CG9231)

Product Code	CSB-CF892906DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VW12
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MLSKTGLIGALVRRSFGTSQMLRETIKNHEPNLERRMLVWTGKYKSQSEIPN FVSQDVM ERCRNKMRIRLANIMIALTAVGCAIMVYSGKQAAKKGESVTKMNLEWHKQFND SQQSEGS APAAK
Source	in vitro E.coli expression system
Target Names	CG9231
Protein Names	Recommended name: UPF0389 protein CG9231
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.