



Recombinant Human Brain protein I3 (BRI3)

Product Code	CSB-CF002811HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95415
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDHKPLLQERPPAYNLEAGQGDIACGPHGYGAIPAAPPPPPYPYLVGTGIPTHH PRVYNIH SRTVTRYPANSIVVVGGCPVCRVGVLEDCFTFLGIFLAILFPFGFICCFALRKRR CPNC GATFA
Source	in vitro E.coli expression system
Target Names	BRI3
Protein Names	Recommended name: Brain protein I3 Alternative name(s): pRGR2
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.