



Recombinant Escherichia coli Protein trbF (trbF)

Product Code	CSB-CF322747ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15068
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MRENKSNPELKIRSTERDYKYISRITGRYAGLSLVFLTAGIVLWTVMDIIFDACID SWKA DPELNNSSYMWNILYAIPTYLYALAAGFLVTFSSVPNVRINIRKYRDIPAEMSYA PGEH IKGGQE
Source	in vitro E.coli expression system
Target Names	trbF
Protein Names	Recommended name: Protein trbF
Expression Region	1-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.