



Recombinant Pan troglodytes NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 (NDUFB6)

Product Code	CSB-CF605274EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQE1
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	TGYTPDEKLRLQQLRELRRRWLKDQELSPREPVLPQKMGPMEKFWNKFLE NKSPWRKMV HGVYKKSIFVFTHVLVPVWIIHYYMKYHVSEKPYGIVEKKSRIFFPGDTILETGEVI PPMK EFPDQHH
Source	in vitro E.coli expression system
Target Names	NDUFB6
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 Alternative name(s): Complex I-B17 Short name= CI-B17 NADH-ubiquinone oxidoreductase B17 subunit
Expression Region	2-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.