



Recombinant Arabidopsis thaliana PRA1 family protein C (PRA1C)

Product Code	CSB-CF627345DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1G3K7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MIFRTNYIVIFIVSIFISMLWQPVHLSVVFVILIVAWLYVYSRDNEPWVIFGVIDDS TLV LVLLVLTIGIFLLTDVSRGIVIGVLAGLPVVLVHGMCRRNTEMLFVLEDDEEKVA MNTSS SSLSSSS
Source	in vitro E.coli expression system
Target Names	PRA1C
Protein Names	Recommended name: PRA1 family protein C Short name= AtPRA1.C
Expression Region	1-127
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.