



Recombinant Mouse NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 (Ndufb6)

Product Code	CSB-CF662371MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3UIU2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	SGYTPDEKLRLQLRELRRRWLKDQELSPREPVLPPRRMWPLERFWDNFLR DGAVWKNMV FKAYRSSLFAVSHVLIPMWFVHYYVKYHMATKPYTIVSSKPRIFPGDTILETGEV IPPMR DFPDQHH
Source	in vitro E.coli expression system
Target Names	Ndufb6
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6 Alternative name(s): Complex I-B17 Short name= CI-B17 NADH-ubiquinone oxidoreductase B17 subunit
Expression Region	2-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.