



Recombinant Danio rerio Succinate dehydrogenase [ubiquinone] cytochrome b small subunit B, mitochondrial (sdhdb)

Product Code	CSB-CF734666DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q68FN7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	AVQQKDHDCSYLISARIHATPSNYAGSGSKAATMHWTGERILSIALLSLAPVAY FCPSA VDYSLAAALTLHGHWGLGQVVTDYVHGDAKIKMANAGLFLSTVTFAGLCYFN YHVDGIC KAVALLWSK
Source	in vitro E.coli expression system
Target Names	sdhdb
Protein Names	Recommended name: Succinate dehydrogenase [ubiquinone] cytochrome b small subunit B, mitochondrial Short name= CybS-B Alternative name(s): Succinate dehydrogenase complex subunit D-B Succinate-ubiquinone oxidoreductase cytochrome b smal
Expression Region	30-158
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.