



Recombinant *Mortierella alpina* Cytochrome b5

Product Code	CSB-CF897598MUM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Y706
Product Type	Transmembrane Protein
Immunogen Species	<i>Mortierella alpina</i> (Oleaginous fungus) (<i>Mortierella renispora</i>)
Sequence	MAELKSFTLADLSQHTTKDSLAIHGKIVYDCTGFIDEHPGGEEVLIDEAGRDA TESFED VGHSDEARDIMSKLLVGEFKTDSSEKPKAKSPSSSTPRPIAAEPSDSGSLQY VLALAVV AGCVIWKVLL
Source	in vitro E.coli expression system
Protein Names	Recommended name: Cytochrome b5
Expression Region	1-130
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.