



Recombinant Escherichia coli O127:H6 Fumarate reductase subunit C (frdC)

Product Code	CSB-CF485659EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7UPX4
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Sequence	MTTKRKPYVRPMTSTWWKKLPFYRFYMLREGTAVPAVWFSIELIFGLFALKNG PEAWAGF VDFLQNPVIVIINLITLAAALLHTKTFELAPKAANIIVKDEKMGPEPIIKSLWAVT VVA TIVILFVALYW
Source	in vitro E.coli expression system
Target Names	frdC
Protein Names	Recommended name: Fumarate reductase subunit C Alternative name(s): Fumarate reductase 15 kDa hydrophobic protein
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.