



Recombinant Human Interferon-induced transmembrane protein 2 (IFITM2)

Product Code	CSB-CF011025HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q01629
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNHIVQTFSPVNSGQPPNYEMLKEEQEVAMLGVPHNPAPPMSTVIHIRSETSV PDHVVWS LFNTLFMNTCCLGFIAFAYSVKSRDRKMVGDTVGAQAYASTAKCLNIWALILGI FMTILL IIPVLVVQAQR
Source	in vitro E.coli expression system
Target Names	IFITM2
Protein Names	Recommended name: Interferon-induced transmembrane protein 2 Alternative name(s): Interferon-inducible protein 1-8D
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.