



Recombinant UPF0041 protein F53F10.3 (F53F10.3)

Product Code	CSB-CF520761CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O01578
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MGPHRLLYQALCKAGDKFVYPVLPFAFAKPAWNHAAGPKTVFFWAPTIKWTLIG AGLADLA RPADKLSLYQNSALFATGAIWTRYCLVITPINYYLSSVNVFFVMCTGLAQLCRIA HYRYQN PDWETKEIMETHN
Source	in vitro E.coli expression system
Target Names	mpc-2
Protein Names	Recommended name: UPF0041 protein F53F10.3
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.