



Recombinant *Synechococcus elongatus* NAD(P)H-quinone oxidoreductase subunit 3 (ndhC)

Product Code	CSB-CF813125FPY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8KPP6
Product Type	Transmembrane Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Sequence	MAKWNSDCSLEFGVFLNGYEYLLGFLLISSLPILSLTASRLLRPGRRGPERR TTYESG MEPIGGAWIQFNVRYYMFALVFVIFDVETVFLYPWAVAFNRLGLLAFVEALIFITI LVVG LAYAWRKGALEWS
Source	in vitro E.coli expression system
Target Names	ndhC
Protein Names	Recommended name: NAD(P)H-quinone oxidoreductase subunit 3 EC= 1.6.5.- Alternative name(s): NAD(P)H dehydrogenase subunit 3 NADH-plastoquinone oxidoreductase subunit 3 NDH-1 subunit 3 Short name= NDH-C
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.