



Recombinant Chicken NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 4 (NDUFB4)

Product Code	CSB-CF015651CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48306
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MASGPSRTAAEEYRPNRYVSLPAELDPATYDTPLEKRRAEAERLAIRARLKRQ YLLQLNT PKPPRVIEDPALLRWDYARTHNVYPTFRPTPKTSFLGAVFAIGPILFWIAAFKTE RVSSA GGVEHQGAPGCKT
Source	in vitro E.coli expression system
Target Names	NDUFB4
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 4 Alternative name(s): Complex I-B15 Short name= CI-B15 GGHPW NADH-ubiquinone oxidoreductase B15 subunit Protein Walter
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.