



Recombinant *Saccharomyces cerevisiae* Mitochondrial outer membrane protein OM14 (OM14)

Product Code	CSB-CF411326STA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6ZLG8
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Sequence	MSATAKHDSNASPNSDSE ₁ DGHHNNKKECAIEYLKARLNSASAVACGYLQAF VSKTQDFA KVCFLELQNPVVLV ₁ NLLLHSSVVCYLCNGYANHNARFLKGGKPNSTVLATTAGA LGLLTLD GIISKKYYSRYDKK
Source	in vitro <i>E.coli</i> expression system
Target Names	OM14
Protein Names	Recommended name: Mitochondrial outer membrane protein OM14 Alternative name(s): Outer membrane protein of 14 kDa
Expression Region	1-134
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.