



Recombinant Arabidopsis thaliana Uncharacterized protein At5g65660 (At5g65660)

Product Code	CSB-CF862870DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9LSK9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MENTQDFSPPHMDASRPSLGFPLGTALLLIIFSLSGIFSCCYHWDKHRSLRRS LANGRP SADIESNPYKPKPPFPEMKKQNLVPLMPGDNTPKFIALPCPCAPPRPEKLT VDVQTP PQSPPVKPARFPVPLY
Source	in vitro E.coli expression system
Target Names	At5g65660
Protein Names	Recommended name: Uncharacterized protein At5g65660
Expression Region	1-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.