



Recombinant *Saccharomyces cerevisiae* ER-derived vesicles protein ERV14 (ERV14)

Product Code	CSB-CF345265SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53173
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	GAWLFILAVVNCINLFGQVHFTILYADLEADYINPIELCSKVNKLITPEAALHGA LSSL FLLNGYWVFLNLPVLAYNLNKIYNKVQLLDATEIFRTLKHKRESFLKLGFHL LMFFF YLYRMIMALIAESGDDF
Source	in vitro E.coli expression system
Target Names	ERV14
Protein Names	Recommended name: ER-derived vesicles protein ERV14
Expression Region	2-138
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.