



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 43, mitochondrial (AIM43)

Product Code	CSB-CF480954SVH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5VT63
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Sequence	IKSLEDLANLDSLDGVDTTELIRDLINEHTTKLNIKKELDMLKKFSQEEESGHEIPV KRFI RPLWMFILMGSSVYLLLFHFSWWKLEHEERESQLKKEVEILEHQLNELIVQDKT HNTSRGK GSNESTHMKPWYRRWFW
Source	in vitro <i>E.coli</i> expression system
Target Names	INA17
Protein Names	Recommended name: Altered inheritance of mitochondria protein 43, mitochondrial Alternative name(s): Found in mitochondrial proteome protein 14
Expression Region	46-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.