



Recombinant Gorilla gorilla gorilla NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 11 (NDUFA11)

Product Code	CSB-CF603880GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQB9
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	APKVFRQYWDIPDGTDCCHRKAYSTTSIASVAGLTAAAYRVTLNPPGTFLEGVA KVGQYTF TAAAVGAVFGLTTCISAHVREKPDDPLNYFLGGCAGGLTLGARTHNYGIGAAA CVYFGIA ASLVKMGQLEGWEVFAKPKV
Source	in vitro E.coli expression system
Target Names	NDUFA11
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 11 Alternative name(s): Complex I-B14.7 Short name= CI-B14.7 NADH-ubiquinone oxidoreductase subunit B14.7
Expression Region	2-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.