



Recombinant Pig Succinate dehydrogenase cytochrome b560 subunit, mitochondrial (SDHC)

Product Code	CSB-CF020908PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D0VWV4
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	LGTTAKEEMERFWNKNLGSNRPLSPHITIYRWSLPMAMSICHRGTGIALSAGV SLFGLSA LLLPGNFESHLELVKSLCLGPTLIYTAKFGIVFPLMYHTWNGIRHLIWDLGKGLTI PQLT QSGVVVLILTVLSSVGLAAM
Source	in vitro E.coli expression system
Target Names	SDHC
Protein Names	Recommended name: Succinate dehydrogenase cytochrome b560 subunit, mitochondrial Alternative name(s): Succinate-ubiquinone oxidoreductase cytochrome B large subunit Short name= CYBL
Expression Region	30-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.