



Recombinant Enterobacteria phage P2 Protein lysA (lysA)

Product Code	CSB-CF346710EDE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51769
Product Type	Transmembrane Protein
Immunogen Species	Escherichia phage P2 (Bacteriophage P2)
Sequence	MKKLSLSLMLNVSLALMLALS LIYPQSVAVNFVAAWAILATVICVVAGG VGVYA TEYVLE RYGRELPPE SLAVKIVTSLFLQPVPWRRRAAALVVVVATFISLVAAGWIFTALIY LVVSL FFRLIRKACRQRLEGREPCQG
Source	in vitro E.coli expression system
Target Names	lysA
Protein Names	Recommended name: Protein lysA
Expression Region	1-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.