



Recombinant Danio rerio Protein C19orf12 homolog (si:ch211-238a12.2, zgc:112052)

Product Code	CSB-CF628658DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1L990
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MPPHVDDVMKLCCELSANQQVKTAVKQSGKGAAAAGGLAFAGGLIGGPLGIA VGGAVGGL LGCWMKSGQFKPLPQVIMELTPDQQARLYEDIVAILGSITWTDVAQLTALVMG NASLQQQ VTAALLSYIHKELQAEVHYID
Source	in vitro E.coli expression system
Target Names	si:ch211-238a12.2
Protein Names	Recommended name: Protein C19orf12 homolog
Expression Region	1-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.