



Recombinant Danio rerio Transmembrane protein 242 (tmem242)

Product Code	CSB-CF690155DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4V8S3
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSVEQASASGAVSDIVEDDKSHLIKGGAFATVATAGMIAGFGATLAVAKKSP DWFNKG IIGSAAVPESGASLALRALGWGSLYAWCGVGLLSLTIWKAMGVHSLQEFRQKM QSIFPAI PKNEDAQANSVPFDWNSIFKSK
Source	in vitro E.coli expression system
Target Names	tmem242
Protein Names	Recommended name: Transmembrane protein 242
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.