



Recombinant Pongo pygmaeus NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13 (NDUFA13)

Product Code	CSB-CF608074EXP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQ88
Product Type	Transmembrane Protein
Immunogen Species	Pongo pygmaeus (Bornean orangutan)
Sequence	GASKVKQDMPPPGGYGPIDYKRNLPRRGLSGYSMLALGIGTLIYGHWSMMK WNRERRRLQ IEDFEARIALPLLQAETDRRTLQMLRENLEEEAIIMKDVPDWKVGESVFHTTR WVAPLI GELYGLRTTEEALHASHGFMWYA
Source	in vitro E.coli expression system
Target Names	NDUFA13
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 13 Alternative name(s): Complex I-B16.6 Short name= CI-B16.6 NADH-ubiquinone oxidoreductase B16.6 subunit
Expression Region	2-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.