



Recombinant *Ashbya gossypii* Peroxisome assembly protein 22 (PEX22)

Product Code	CSB-CF765636DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q758W8
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	MVSKASRDQLRKYGAVSLASLLVAASIVAYRWWNAAPSIEVEKLLRRSVSRCV VVTQGIQ NEDMIHDLLFEDTVMLLAPGCTAEGRLKSASRENAYKVISCTTWQSVWACVR HFRKHTLL VRTSEVPSGVPADIGGYVSDISDI
Source	in vitro <i>E.coli</i> expression system
Target Names	PEX22
Protein Names	Recommended name: Peroxisome assembly protein 22 Alternative name(s): Peroxin-22
Expression Region	1-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.