



Recombinant Arabidopsis thaliana Probable gamma-secretase subunit PEN-2 (At5g09310)

Product Code	CSB-CF872341DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FY84
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MEATRSDDP SLNPIRNRNPNPNPNPLSTI ISSAQVWPTIDG PLGLTEEASVD YARRFY KFGFALLPWLW FVNCFYFWPVL RHSRAFPQIRNYV VRSAIGFSVFTALLSAWA LTFSIGG EQLFGPLYDKLVMY NVADRLGLSGLA
Source	in vitro E.coli expression system
Target Names	At5g09310
Protein Names	Recommended name: Probable gamma-secretase subunit PEN-2
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.