



Recombinant Arabidopsis thaliana Protein PLANT CADMIUM RESISTANCE 9 (PCR9)

Product Code	CSB-CF317388DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CW98
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSEQEGKNEKKVTEGQWTTGLYDCLSEDISTCCFTWVCPCVAFGRGIAEILDKG ETSRGLA GLMVMAMSSIGCGWYYASKYRAKLRHQYALPEAPCADGAIHCFCCPCALTQE HRELKHRG LDPSLGWNIENGLNSNTPPFVASGMDR
Source	in vitro E.coli expression system
Target Names	PCR9
Protein Names	Recommended name: Protein PLANT CADMIUM RESISTANCE 9 Short name= AtPCR9
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.