



Recombinant Human Oligosaccharyltransferase complex subunit OSTC (OSTC)

Product Code	CSB-CF882089HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NRP0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METLYRVPFLVLECPNLKLLKPPWLHMPSAMTVYALVVVSFYFLITGGIYDVIVE PPSVG SMTDEHGHQRPVAFLAYRVNGQYIMEGLASSFLFTMGGLGFILDRSNAPNIPK LNRFLLLFIGFVCVLLSFFMARVFMRMKLPGYLMG
Source	in vitro E.coli expression system
Target Names	OSTC
Protein Names	Recommended name: Oligosaccharyltransferase complex subunit OSTC Alternative name(s): Hydrophobic protein HSF-28
Expression Region	1-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.