



Recombinant Mouse Amphiregulin (Areg)

Product Code	CSB-CF001986MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31955
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	VIKPKKNKTEGEKSTEKPKRKKKGGKNGKGRRNKKNPCTAKFQNF CIHGE CRYIENLE VVTCNCHQDYFGERCGEKSMKTHSEDDKDL SKIAVVAVTIFVSAILLAAIGIGIVI TVHL WKRYFREYEGETEERRRLRQENGTVHAIA
Source	in vitro E.coli expression system
Target Names	Areg
Protein Names	Recommended name: Amphiregulin Short name= AR Alternative name(s): Schwannoma-derived growth factor Short name= SDGF
Expression Region	100-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.