



Recombinant Escherichia coli O17:K52:H18 UPF0208 membrane protein YfbV (yfbV)

Product Code	CSB-CF484557EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7N5Q7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Sequence	MSTPDNRSVNFFSLFRRGQHYSKTWPLEKRLAPVFVENRVIKMTRYAIRFMPP IAVFTLC WQIALGGQLGPAVATALFALS LPMQGLWWLGKRSVTPLPPAILNWFYEV RGK LQESGQVL APVEGKPDYQALADTLKRAFKQLDKTFLDDL
Source	in vitro E.coli expression system
Target Names	yfbV
Protein Names	Recommended name: UPF0208 membrane protein YfbV
Expression Region	1-151
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.