



Recombinant Escherichia coli O8 UPF0266 membrane protein yobD (yobD)

Product Code	CSB-CF485978E00
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7M294
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Sequence	MTITDLVLILFIAALLAFAIYDQFIMPRRNGPTLLAIPLLRGRIDSVIFVGLIVILIYN NVTNHGALITTWLLSALALMGFYIFWIRVPKIIFKQKGGFFANVWIEYSRIKAMNL SEDG VLVMQLEQRLLIRVRNIDDLEKIYKLLVSTQ
Source	in vitro E.coli expression system
Target Names	yobD
Protein Names	Recommended name: UPF0266 membrane protein yobD
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.