



Recombinant Escherichia coli O127:H6 UPF0266 membrane protein yobD (yobD)

Product Code	CSB-CF487357EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7USJ6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Sequence	MTITDLVLILFIAVLLAFAIYDQFIMPRRNGPTLLAIPLLRGRIDSVIFVGLIVILIYN SVTNHGALITTWLLSALALMGFYIFWIRVPKIIFKQKGGFFANVWIEYSRIKAMNL SEDG VLVMQLEQRLLIRVRNIDDLEKVYKLLVSTQ
Source	in vitro E.coli expression system
Target Names	yobD
Protein Names	Recommended name: UPF0266 membrane protein yobD
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.