



Recombinant Human Chemokine-like factor (CKLF)

Product Code	CSB-CF005458HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UBR5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDNVQPKIKHRPFCFSVKGHVKMLRLALTVTSMFFIIAQAPEPYIVITGFEVTVI LFFI LLYVLRDLRLMKWLFWPLLDIINSLVTTVFMLIVSVLALIPETTTTLTVGGGVFALV TAVC CLADGALIYRKLLFNPSGPYQKKPVHEKKEVL
Source	in vitro E.coli expression system
Target Names	CKLF
Protein Names	Recommended name: Chemokine-like factor Alternative name(s): C32
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.