



Recombinant Flagellar protein FliL (fliL)

Product Code	CSB-CF364716EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ABX9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MTDYAISKKSKRSLWIPILVFITLAACASAGYSYWHSHQVAADDKAQQRVVPSP VFYALD TFTVNLGDADRVLYIGITLRLKDEATRSRLSEYLPEVRSRLLLLFSRQDAAVLAT EEGKK NLIAEIKTTLSTPLVAGQPKQDVTDVLYTAFILR
Source	in vitro E.coli expression system
Target Names	fliL
Protein Names	Recommended name: Flagellar protein FliL
Expression Region	1-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.