



Recombinant Human RING finger protein 122 (RNF122)

Product Code	CSB-CF019825HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H9V4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MHPFQWCNGCFCGLGLVSTNKSCSMPPISFQDLPLNIYMVIFGTGIFVFMLSLI FCCYFI SKLRNQAQSERYGYKEVVVLKGDAAKLLQLYGQTCAVCLEDFKGGKDELGVLPCQ HAFHRKCL VKWLEVRCVCPMCNKPIASPSEATQNIGILLDELV
Source	in vitro E.coli expression system
Target Names	RNF122
Protein Names	Recommended name: RING finger protein 122
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.