



Recombinant Human Uncharacterized protein C10orf111 (C10orf111)

Product Code	CSB-CF843145HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N326
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MESLQTPQHRENQDKREKEYGVKHMPMGNNAGNLEPEKRKAVRVALSSATA AQNIPSSVH CGCSKQWRLRLPSESLQSRGQVMKRPNNILKLRNLDLLIYPWPELRRRQVAS DLMSLLLL PAFSGLTWAPFLFLFTYLPPFLNLLTVGFVSYFLV
Source	in vitro E.coli expression system
Target Names	C10orf111
Protein Names	Recommended name: Uncharacterized protein C10orf111
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.