



Recombinant Mouse Protein FAM162A (Fam162a)

Product Code	CSB-CF880508MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9D6U8
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MWSLGGLRLAAGHCLRLYERNASSSLRFTRNTDLKRINGFCTKPKQESPKTPT QSYRHGVP LHKPTDFEKKILLWSGRFKKEEEIPETISFEMLDAAKNKLRVKVSYLMIALTVAG CIYMV IEGKKAAGRHESLTSNLNLERKARLREEAAMKAKTD
Source	in vitro E.coli expression system
Target Names	Fam162a
Protein Names	Recommended name: Protein FAM162A Alternative name(s): E2-induced gene 5 protein homolog
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.