



Recombinant Zea mays Oleosin Zm-I (OLE16)

Product Code	CSB-CF321111ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13436
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	ADHHRGATGGGGGYGDLQRGGGMHGEAQQQQKQGAMMTALKAATAATFG GSMLVLSGLIL AGTVIALTVATPVLVIFSPVLVPAIALALMAAGFVTSGGLGVAALSVFSWYKY LTGKH PPAADQLDHAKARLASKARDVKDAAQHRIDQAQGS
Source	in vitro E.coli expression system
Target Names	OLE16
Protein Names	Recommended name: Oleosin Zm-I Alternative name(s): Lipid body-associated major protein Lipid body-associated protein L3 Oleosin 16 kDa
Expression Region	2-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.