



Recombinant Rat Melatonin receptor type 1A (Mtnr1a)

Product Code	CSB-CF015189RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49218
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	CYICHSLKYDRIYSNKNSLCYVFLIWLTLIAIMPNLQTGTLQYDPRIYSCTFTQS VSSA YTIALVVFHFVVPMIIVTFCYLRIWILVLQVRRRVKPDSPKPKLKPQDFRNFVTMF VVFVL FALCWAPLNFIGLIVASDPATMAPRIPEWLFVASY
Source	in vitro E.coli expression system
Target Names	Mtnr1a
Protein Names	Recommended name: Melatonin receptor type 1A Short name= Mel-1A-R Short name= Mel1a receptor
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.