



Recombinant Tetraodon nigroviridis FUN14 domain-containing protein 1 (fundc1)

Product Code	CSB-CF700192TCA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4RY26
Product Type	Transmembrane Protein
Immunogen Species	Tetraodon nigroviridis (Spotted green pufferfish) (Chelonodon nigroviridis)
Sequence	MATVDHREAGQDDPESEDEVYEVVDLLEYARRHQWWSRVFGSNSGPIAEKY SVTTQLVMG GLTGWCAGYLFQRVVGKIAATAVGGGFLLLQIANHSGYVQVDWKKVEKDVNKA KKRLKKKA NKA VPEINTLIEEATDFIKKNIVLSGGFVGGFLLGLAS
Source	in vitro E.coli expression system
Target Names	fundc1
Protein Names	Recommended name: FUN14 domain-containing protein 1
Expression Region	1-158
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.